

THE CYCLE OF COLOR


Making sure things don't get too colorful when using plastic recyclates



MEASUREMENT OF THE COLOR OF PLASTIC RECYCLATES BY MEANS OF THE LABORATORY DEVICE SPECTRO-3-0°/45°-MSM-ANA-LAB-P



Color measurement of plastic recyclate samples with the lab colorimeter SPECTRO-3-0°/45°-MSM-LAB-P



Lab colorimeter SPECTRO-3-0°/45°-MSM-LAB-ANA-P with cab eos2 label printer and RAL - plastic cards as well as recyclate samples




Evaluation unit including sensors, panel PC, USB interfaces, calibration card holder, support for sample holder, pellet hopper with sight glass, slider



Calibration card holder, holder for pellet samples, pellet hopper with sight glass, slider

Label printer



Calibration cards




Pellet samples




Pellet funnel with sight glass


Laboratory color measurement system SPECTRO-3-0°/45°-MSM-ANA-P



Electrical and mechanical interfaces of the lab colorimeter SPECTRO-3-0°/45°-MSM-LAB-ANA-P




- glass plate
- calibration card holder
- holder
- calibration card
- mechanical interface for pellet hopper
- slider
- holder for pellet sample container






- pellet hopper with sight glass
- recyclate
- pellet sample container


Calibration card holder with inserted RAL - plastic card (left) as well as with flanged pellet hopper and pellet sample container (right)






RAL - plastic cards with indication of L*a*b* - values, measured according to the diffuse/8° - and 45°/0° - method. Up to 300 different RAL - plastic cards are available.



For the calibration of the lab color measuring system the RAL - plastic cards are used. Three different white tones are available for the white balance: RAL9003-P, RAL9010-P and RAL9016-P. One of these three RAL - card types should be used for the white balance. After placing one of the three cards in the calibration card holder, the software of the measuring system adjusts to the white surface of the RAL - plastic card in an optimal way.



In the following example the RAL9003-P was used as reference card for the white balance. All L*a*b* - color values of the approx. 300 RAL - plastic cards are available on the files **RAL Plastic Plate measured with Hunter Lab - 45°_0°** as well as **RAL Plastic Plate Standard d_8°**.



Selection of the color measurement standard: $45^{\circ}/0^{\circ}$ or $d/8^{\circ}$ method; $45^{\circ}/0^{\circ}$ is recommended here, since the laboratory color measurement system also works according to this measurement method. After selecting the file the used RAL - card type must be entered, in the example here: RAL9003-P.

SPECTRO3 MSM DOCAL Scope V1.0

SPECTRO3 MSM DOCAL Scope V1.0

CONNECT | PARA1 | TEACH | CALIB | DOCU

START CALIBRATION ASSISTANT

REFERENCE

COLOR SPACE OF REF SETVALUE L* a* b*

SETVALUE TABLE

C SPACE L* a* b* ROWS 1

GET CF

RAM SEND GO

EE GET STOP

FILE

CALIB REF | CALIB DATA | SENSOR DATA

Scroll all CLOSE

SETVALUE TABLE

	L*	a*	b*
1	91.9100	-0.2700	0.5200
2	0.0000	0.0000	0.0000
3	0.0000	0.0000	0.0000
4	0.0000	0.0000	0.0000
5	0.0000	0.0000	0.0000
6	0.0000	0.0000	0.0000
7	0.0000	0.0000	0.0000
8	0.0000	0.0000	0.0000
9	0.0000	0.0000	0.0000
10	0.0000	0.0000	0.0000
11	0.0000	0.0000	0.0000
12	0.0000	0.0000	0.0000
13	0.0000	0.0000	0.0000
14	0.0000	0.0000	0.0000
15	0.0000	0.0000	0.0000
16	0.0000	0.0000	0.0000
17	0.0000	0.0000	0.0000
18	0.0000	0.0000	0.0000
19	0.0000	0.0000	0.0000
20	0.0000	0.0000	0.0000
21	0.0000	0.0000	0.0000
22	0.0000	0.0000	0.0000
23	0.0000	0.0000	0.0000
24	0.0000	0.0000	0.0000
25	0.0000	0.0000	0.0000

ACTUAL XYZ VALUE TABLE

	X SI	Y SI	Z SI
1	2761	2956	3199
2	0	0	0
3	0	0	0
4	0	0	0
5	0	0	0
6	0	0	0
7	0	0	0
8	0	0	0
9	0	0	0
10	0	0	0
11	0	0	0
12	0	0	0
13	0	0	0
14	0	0	0
15	0	0	0
16	0	0	0
17	0	0	0
18	0	0	0
19	0	0	0
20	0	0	0
21	0	0	0
22	0	0	0
23	0	0	0
24	0	0	0
25	0	0	0

RAW X 2761

RAW Y 2956

RAW Z 3199







Message!

Do you want to add another card?






Yes No

18:46
28.12.2021


1st calibration card



After calibration to the white reference (RAL9003-P), several local reference points can now be stored in the color space for more precise adjustment of the color measurement system. It is recommended to use RAL - plastic cards, which visually corresponds to the produced recycle in terms of color. After selecting a suitable RAL plastic card, it is entered into the software (here: RAL9002-P) and then positioned in the calibration card holder. After clicking the OK software button twice, the $L^*a^*b^*$ color values of the RAL plastic card and the XYZ color values of the color measuring system are entered into the respective calibration table.



The calibration table can be supplemented by entering further reference points (here: RAL7035-P). For this purpose, the RAL plastic card RAL7035-P is placed in the calibration card holder. Enter "RAL7035-P into the software and then confirm twice with a mouse click on the OK software button.






Adding another calibration card (here: RAL5024-P) to the calibration table:

*Place the RAL - card in the calibration card holder.

*Input of "RAL5024-P" into the software.

*Confirm twice by mouse click with OK.



Again the software inquires for a further RAL plastic card, in our example this is the RAL9011-P. After inserting the RAL plastic card into the calibration card holder and entering "RAL9011-P" in the software and confirming twice by clicking the OK software button, an entry is made, in our example the 5th in the SETVALUE TABLE as well as in the ACTUAL XYZ VALUE TABLE.

SPECTRO3 MSM DOCAL Scope V1.0

CONNECT PARA1 TEACH CALIB DOCU

START CALIBRATION ASSISTANT

REFERENCE
COLOR SPACE OF REF SETVALUE L*a*b*

SETVALUE TABLE
C SPACE L*a*b* ROWS 5

GET CF

RAM SEND GO
EE GET STOP
FILE

CALIB REF CALIB DATA SENSOR DATA

Scroll all CLOSE

SETVALUE TABLE

	L*	a*	b*
1	91.9100	-0.2700	0.5200
2	84.5900	-1.0900	5.5400
3	78.6600	-1.0100	0.9600
4	56.2800	-7.5700	-19.6900
5	14.2500	0.0500	-0.9500
6	0.0000	0.0000	0.0000
13	0.0000	0.0000	0.0000
14	0.0000	0.0000	0.0000
15	0.0000	0.0000	0.0000
16	0.0000	0.0000	0.0000
17	0.0000	0.0000	0.0000
18	0.0000	0.0000	0.0000
19	0.0000	0.0000	0.0000
20	0.0000	0.0000	0.0000
21	0.0000	0.0000	0.0000
22	0.0000	0.0000	0.0000
23	0.0000	0.0000	0.0000
24	0.0000	0.0000	0.0000
25	0.0000	0.0000	0.0000

ACTUAL XYZ VALUE TABLE





	X SI	Y SI	Z SI
1	2761	2956	3199
2	2210	2381	2379
3	1894	2038	2183
4	770	920	1492
5	61	66	79
6	0	0	0
7	0	0	0
8	0	0	0
9	0	0	0
10	0	0	0
11	0	0	0
12	0	0	0
13	0	0	0
14	0	0	0
15	0	0	0
16	0	0	0
17	0	0	0
18	0	0	0
19	0	0	0
20	0	0	0
21	0	0	0
22	0	0	0
23	0	0	0
24	0	0	0
25	0	0	0

RAW X 61
RAW Y 66
RAW Z 79


Message!
Do you want to add another card?
Yes No

18:52
28.12.2021


5th calibration card




In our example we are finished with the RAL plastic cards to which the color measuring system should be calibrated and the question for further RAL plastic cards can be answered with "No" (mouse click on the "No" - field of the software). The software informs that the calibration was successful and that the calibration data has been taken over and stored in the non-volatile memory of the controller in the color measuring device after confirming "OK" twice via mouse click. Furthermore, the software informs that the CALIB mode UCAL has been activated, i.e. the sensor system now uses the currently entered calibration values to determine the color values. The activation or the UCAL mode was stored in the non-volatile memory of the sensor-internal controller.



PARA1 user interface of the SPECTRO3 MSM DOCAL Scope V1.0 software after calibration. The value in the AVERAGE input field should be changed from 1 to 32768. Then click on EE and confirm with SEND via mouse click.



The actual measurement process is now to be carried out on the basis of a recyclate that is very close in color to the RAL plastic card RAL5024-P from a purely visual point of view. In this process, several recyclate samples are taken at specific time intervals from ongoing production using the sample containers and fed to the lab colorimeter. First, it is necessary to attach or flange the recyclate pick-up unit (pellet hopper unit) to the calibration card pick-up unit. Four elongated holes are provided on the rear of the recyclate pick-up unit for this purpose. In addition, four bars matching the elongated holes are provided on the calibration card pick-up unit to allow the recyclate pick-up unit to be flanged on. After the recyclate pick-up unit has been mounted, the recyclate sample can be fed in. After filling the recyclate pick-up unit, the empty recyclate container should be positioned below the recyclate pick-up unit (pellet hopper unit) in the recess provided for this purpose.




	a*	b*	L*	deltaE				
0	0.00	0.00	0.00	0.00	0.00	0.00	0.00	Red
1	0.00	0.00	0.00	0.00	0.00	0.00	0.00	Green
2	0.00	0.00	0.00	0.00	0.00	0.00	0.00	Blue

TEACH DATA TO No.: 0 Inc [checkbox] [Power icon]






	a*	b*	L*	deltaE				
0	-0.25	6.47	85.69	10.00	0.00	0.00	0.00	Red
1	-7.22	-20.24	61.69	10.00	0.00	0.00	0.00	Green
2	0.00	0.00	0.00	0.00	0.00	0.00	0.00	Blue

The measurement can now begin. So that the deviations dL^* , da^* and db^* from a reference value can also be displayed, this must be first defined. This is made possible by entering the $L^*a^*b^*$ reference color values in the so-called TEACH TABLE. The SPECTRO-3-0°/45°-MSM-LAB-ANA-P software can store up to three different reference color triplets. The values can be entered via the keyboard after clicking on the respective field in the TEACH TABLE, but also by means of suitable recycle samples, whose values (after the recycle sample has been filled into the recycle pick-up unit, named as the pellet hopper unit) are stored in the $a^*b^*L^*$ fields of the selected line after selecting the line number (in the TEACH TABLE) and after clicking on the TEACH DATA TO field of the software.




The TEACH TABLE is now filled in and the color measurement of the recycle samples (in our example the pastel blue recycle samples) can be continued. The pastel blue recycle reference was stored in row 1 of the TEACH TABLE: $a^* = -7.22$, $b^* = -20.44$, $L^* = 61.69$. After clicking the GO software button, the current $L^*a^*b^*$ values of the recycle sample, as well as the $dL^*da^*db^*$ color deviation values (deviation to the reference values), and furthermore the dE value, thus the deviation from the reference, stored in the TEACH TABLE, in the color space in entirety, are displayed in the column roughly in the middle of the PC screen.



The color values $L^*a^*b^*$ of the respective recyclate sample, as well as their deviations $dL^*da^*db^*$ from the matching reference in the TEACH TABLE, can now be saved in a file that has yet to be created. After clicking the software button "SELECT RECORD FILE" the file name and the file location can be defined, in our example: \\DOCAL\RecordFile_PASTEL_BLUE_RECYCLATE_038_II.dat. After placing the recyclate sample in the recyclate pick-up unit (pellet hopper unit), the measurement can now be triggered by clicking the "CAPTURE DATA FRAME" software button. Parallel to the recording of the color values $L^*a^*b^*$, the color deviations $dL^*da^*db^*$ to the matching color reference $a^*b^*L^*$ stored in row1 in the TEACH TABLE, the date as well as the time, a graphical display also takes place in the SPECTRO3 MSM DOCAL Scope V1.0 software, also in the DOCU menu item.



Graphical as well as numerical display of the $L^*a^*b^*$ - color values, the $dL^*da^*db^*$ - color deviations and the dE - value for the matching color reference in the TEACH TABLE in the menu item DOCU.



PRINTOUT LINE1 (Max 20 characters)
 PASTEL_BLUE_RECYCLATE_010

PRINTOUT LINE2 (Max 20 characters)
 12-25-2021 12:49

FILL LINE2 WITH DATE

PRINT DOCUMENTATION DATA

PRINT DOCUMENTATION DATA

Drucken

Name: CAB EOS2/80 Eigenschaften

Status: Bereit

Typ: CAB EOS2/300

Standort: USB001

Kommentar: Ausgabe in Date

Druckbereich: Alle

Beispiel: Anzahl Exemplare: 1

Seiten von 0 bis 0

Markierung

OK Abbrechen






Message! X


Print successful

OK


A numeric output of the data including text field can be done for example by means of a label printer (cab eos2). For text input, 2 text lines are available (each with a maximum of 20 characters), but the 2nd line can also be used for the date and time (click on "FILL LINE2 WITH DATE"). After filling in the two lines, the label can now be created by clicking on "PRINT DOCUMENTATION DATA" and then selecting the printer, including a mouse click on the OK - software button. The software then acknowledges the successful label printing with "Print successful". With a mouse click on the OK software button, the software and the color measurement system are available for a new measurement.



After the label has been printed, the recyclate sample can be removed from the recyclate pick-up unit (pellet hopper unit) by pulling the slide and at the same time filled into the empty recyclate sample container. The label created for the recyclate sample can then be attached to the sample container.




1. After another pellet sample has been taken, the pellets are filled into the recycle pick-up unit (pellet hopper unit).
2. Placing the now empty pellet sample container in the recess provided for this purpose below the pellet hopper unit.
3. Activation of the measurement by mouse click on the CAPTURE DATA FRAME software button.
4. In the right part of the software interface the current $L^*a^*b^*$ - values and the $dL^*da^*db^*$ - values numerically and the $dL^*da^*db^*$ - values numerically as well as graphically. In addition, the values are saved in the selected file together with the time and date.
5. In line 1 on the label to be created, text but also numbers and special characters can be entered. Line 2 can be used in the same way as line 1, or the date and the time can be entered.
6. By clicking on the OK - software button, the label printer will be prompted to create the label.
7. The label will be created on the label printer.
8. A message is displayed, that the printing process was successful.
9. The pellets can now be removed from the pellet hopper unit by pulling the slide. The pellets are thereby filled into the already placed, empty pellet sample container. Now only the lid has to be screwed on and
10. Peeling the label from the backing material.
11. Sticking the label onto the pellet sample container.



2nd measurement



1.



2.



CAPTURE DATA FRAME

3.



4.



PRINTOUT LINE1 (Max 20 characters)
 PASTEL_BLUE_RECYCLAT_010
 PRINTOUT LINE2 (Max 20 characters)
 12-25-2021 12:49
 FILL LINE2 WITH DATE
 PRINT DOCUMENTATION DATA

5.




6.



7.



8.



9.



10.



11.


3rd measurement



3rd measurement



1.



2.



CAPTURE DATA FRAME

3.



4.



PRINTOUT LINE1 (Max 20 characters)
 PASTEL_BLUE_RECYCLATE_010
 PRINTOUT LINE2 (Max 20 characters)
 12-25-2021 12:49

5.



Drucker
 Name: CAB EOS2-00
 Status: Bereit
 Typ: CAB 5032/300
 Standort: USB001
 Kommentar: Ausgabe in Datei
 Druckbereich:
 Alle
 Markierung
 Exemplare: Anzahl Exemplare: 1

6.




7.



Message X
 Print successful!

8.



9.




10.



11.




100th measurement



100th measurement

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O
1															
2	Record results of: SPECTRO3 MSM DOCAL Scope V1.0														
3															
4	DATE	TIME	X	Y	Z	L*	a*	b*	delta E	delta L*	delta a*	delta b*	COLOR	TEMP	
5	12-25-2021	12:39:10	983	1103	1819	61.688	-7.122	-19.870	0.382	0.000	0.097	0.369		1	18
6	12-25-2021	12:39:16	983	1103	1817	61.688	-7.122	-19.820	0.429	0.000	0.097	0.418		1	18
7	12-25-2021	12:39:18	983	1103	1817	61.688	-7.122	-19.820	0.429	0.000	0.097	0.418		1	18
8	12-25-2021	12:39:20	983	1103	1819	61.688	-7.122	-19.870	0.382	0.000	0.097	0.369		1	18
9	12-25-2021	12:39:21	983	1103	1819	61.688	-7.122	-19.870	0.382	0.000	0.097	0.369		1	18
10	12-25-2021	12:39:22	983	1103	1819	61.688	-7.122	-19.870	0.382	0.000	0.097	0.369		1	18
11	12-25-2021	12:39:23	983	1104	1819	61.717	-7.249	-19.819	0.422	0.029	-0.030	0.420		1	18
12	12-25-2021	12:39:44	986	1106	1823	61.761	-7.084	-19.854	0.414	0.073	0.134	0.385		1	18
13	12-25-2021	12:39:45	986	1106	1823	61.761	-7.084	-19.854	0.414	0.073	0.134	0.385		1	18
14	12-25-2021	12:39:46	986	1106	1823	61.761	-7.084	-19.854	0.414	0.073	0.134	0.385		1	18
15	12-25-2021	12:39:46	986	1106	1823	61.761	-7.084	-19.854	0.414	0.073	0.134	0.385		1	18
16	12-25-2021	12:39:47	986	1106	1823	61.761	-7.084	-19.854	0.414	0.073	0.134	0.385		1	18
17	12-25-2021	12:39:52	987	1108	1824	61.805	-7.178	-19.815	0.442	0.117	0.041	0.424		1	18
18	12-25-2021	12:39:53	987	1108	1824	61.805	-7.178	-19.815	0.442	0.117	0.041	0.424		1	18
19	12-25-2021	12:39:54	988	1109	1827	61.827	-7.144	-19.851	0.419	0.140	0.075	0.388		1	18
20	12-25-2021	12:39:55	988	1109	1827	61.827	-7.144	-19.851	0.419	0.140	0.075	0.388		1	18
21	12-25-2021	12:39:56	987	1108	1825	61.805	-7.178	-19.839	0.418	0.117	0.041	0.399		1	18
22	12-25-2021	12:39:59	987	1109	1828	61.827	-7.273	-19.887	0.382	0.140	-0.054	0.351		1	18
23	12-25-2021	12:40:00	987	1109	1828	61.827	-7.273	-19.887	0.382	0.140	-0.054	0.351		1	18
24	12-25-2021	12:40:00	990	1110	1832	61.857	-7.045	-19.948	0.379	0.169	0.173	0.291		1	18
25	12-25-2021	12:40:01	990	1110	1832	61.857	-7.045	-19.948	0.379	0.169	0.173	0.291		1	18
26	12-25-2021	12:40:02	990	1110	1832	61.857	-7.045	-19.948	0.379	0.169	0.173	0.291		1	18
27	12-25-2021	12:40:03	990	1110	1832	61.857	-7.045	-19.948	0.379	0.169	0.173	0.291		1	18

Open the file RecordFile_PASTEL_BLUE_RECyclAT_038_II



Graphic display of the measured 100 L*a*b* color values